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www.cenbimo.com
www.histosonda.com

Fabricante: CENBIMO S.L.
C/ Doctor Iglesias Otero, s/n
27004 Lugo, Spain

Assay for 1 individual reaction in lyophilized single test vial.

Product classification:

*EU countries: "For In Vitro Diagnostics"
*Non- EU countries: "For Research Use Only"

INTRODUCTION

Alpha Fetal Protein is the most abundant seric protein in the fetus which demonstrates great similarity to the physical properties of serum albumin.

In the adult, the seric concentration of AFP is very low, apart from when the patient presents with a hepatocarcinoma or teratoma.

The genes that code for AFP and serum albumin are syntenic genes that are situated on the long arm of chromosome 4.

After birth the organism inverts the production of AFP so that it declines rapidly and also considerably increases the production of serum albumin.

There are commercial antibodies against AFP for the diagnosis of tumors that produce this protein but they almost always produce so much tissue background that a clear view becomes impossible.

Histosonda Alpha Fetoprotein is a single stranded DNA fragment of 311 nucleotides that is targeted against the RNA of AFP. Its use is centered on the diagnosis of primitive liver tumors and embryonic tumors of the gonads.

Up to 70% of hepatocarcinomas are producers of AFP. The use of this probe combined with Histosonda Serum Albumin permits a correct diagnosis of primitive liver tumors and their differential diagnosis with metastatic tumors in the liver.

The embryonic tumors of the gonads and midline can show an important degree of complexity and distinct levels of differentiation not well recognisable with conventional histological techniques. Histosonda Alpha Fetal Protein permits the identification of the number and localization of cells producing AFP in a tissue section.

INTENDED USE

For use in In Vitro Diagnosis.
Histosonda AFP is useful for the diagnosis of primitive liver tumors and embryonic tumors of the gonads.

WARNINGS AND PRECAUTIONS

Histosonda AFP has been designed for professional use in In Vitro Diagnosis and must be manipulated by qualified and accordingly trained personnel. In order to obtain the best results, the instructions contained in the manual must be followed. Any change to the indicated temperatures, times or any other step of the process can lead to poor results.

COMPONENTS

1 single test tube of lyophilized Histosonda AFP.
Histosonda AFP consists of a segment of single-stranded DNA complementary to expressed RNA with a length of 311 nucleotides. This probe has been labeled with Digoxigenin.

STORAGE CONDITIONS

Supplied reagents are stored at room temperature until expiration date. After being re-

constituted the probe will remain stable at 4°C for two weeks in a DNAase-free environment.
Do not use after expiration date.

SAMPLES

Any paraffin block section in which Alpha Fetal Protein RNA presence is to be studied. Sections of 4-6 micrometers in width are sufficient to conduct the study. Preferably, the cut should be recent (no more than thirty days old). Assay results are not affected by block age. Studies have been carried out in the manufacturers' laboratories using 20 year old paraffin blocks with optimal results.

INTERPRETATION OF RESULTS

Samples in which Alpha Fetal Protein expression is observed will show a brownish color in the cell cytoplasm, which will contrast over the blue-violet background given by hematoxylin staining. The pathologist will evaluate the results according to their experience, drawing conclusions from the staining of the sample, in parallel with the staining observed in the positive and negative controls.

ASSAY LIMITATIONS

Histosonda AFP has been optimized to detect RNA expression in formalin-fixed, paraffin-embedded tissues. Its use is not recommended for other types of samples or preparation techniques.

The correct operation of this product has been validated using the protocols indicated in the instructions manual. The use of other procedures or the modification of the recommended protocols may lead to erroneous results.

The results from this assay must be evaluated by the pathologist in combination with the rest of available patient clinical data. In order to obtain optimal and reproducible results it is important to rigorously maintain the time and temperature conditions indicated in the procedure.

Emission date: 01/03/2010

PROCEDURE

For Digoxigenin- labeled probes

BASIS OF THE METHOD

Chromogenic In Situ Hybridization (CISH) is a technique used to determine the presence of a DNA or RNA sequence or to study gene expression, as well as for the simultaneous evaluation of tissue morphology by white light microscopy. A labeled probe of known sequence hybridizes with its target in the tissue being studied. This hybridization is then detected by an immunohistochemical process.

The Histosondas targets are messenger RNA and they are therefore designed to detect and visualise gene expression.

HistoSonda Protocol

Before starting pre-heat the humid chamber to 62°C.

1. Reconstitute the probe

1. Add 65µl of (good quality) distilled water to the HistoSonda tube, briefly vortex and centrifuge if possible.

2. Deparaffinization

- Heat slides at 62°C for 10mins using a hot plate or incubator.
- Immerse slides in:
 - Xylene: 10 mins
 - Xylene: 5mins
 - Absolute alcohol (ethanol or isopropanol): 1min X 2
 - Alcohol 96%: 1min X 3
 - Methanol containing 0.3% H₂O₂ or 3% H₂O₂ only :5mins
- Wash well with distilled water, leave standing in water.

2b. Inhibition of unspecific DNA binding (optional)

Generally this step is not necessary unless the tissue to be hybridized contains a great quantity of polymorphonuclear eosinophil leukocytes (generally bone marrow and gastric tissues).

- Place slides in 200ml distilled water in a microwave until the liquid boils uniformly (normally around 2mins).
- Immediately transfer slides to distilled water at room temperature.
- Alternatively place slides in boiling distilled water for 30 seconds and immediately transfer to distilled water at room temp.

! Important: excess heat treatment will result in background staining. Remove slides immediately. !

3. Deproteinization

- Reconstitute Proteinase K from Cenbimo with 200µl of PBS (1 vial for each slide).
- Remove excess water from individual slides with tissue paper.
- Cover tissue sections with the 200 µl of Proteinase K and incubate in a humid chamber for exactly 10mins at room temp.
- Wash well with water.
- Transfer to PBS pH 7.4 for 2mins.

! (Important: if the slides have been boiled in the microwave previously only use Proteinase K for 5 minutes as the heat treatment will have significantly sensitized the tissue sections. For bone marrows continue to use 10mins after heat treatment. Bone marrows fixed in formic acid or EDTA require a Proteinase K digestion of 20min at 55°C.) !

4. Incubation with the probe

- Remove excess buffer from sections as described previously.
- Add the 65µl of probe to the tissue section ensuring the entire section is completely covered and avoiding air bubbles.
- Place the slides in a horizontal position in a humid chamber and close well.
- Incubate at 62°C for 1hr.
- If no suitable humid chamber is available drop the 65µl of probe over a 24 X 50 mm coverslip and place over the tissue section. Seal the edges of the coverslip with rubber cement. Incubate in a 62°C incubator or over a hot plate for 1 hour. After incubation carefully remove the rubber cement and coverslip and continue the protocol.

5. Washing the probe

- Wash the surface of the section vigorously with PBS pH 7.4 to remove the probe.
- Agitate in PBS for 5mins.

6. Revealing the probe

Protocol for manual revealing of the probe. Alternatively, an automatic immunohistochemistry apparatus may be used.

- Reconstitute Anti-Digoxin from Cenbimo with 1ml of distilled water. Each vial is valid for 10 slides.
- Remove excess buffer from sections as before.
- Cover sections with 100µl of Anti-Digoxin and incubate in a humid chamber for 30mins at room temp.
- Wash vigorously with PBS, pH 7.4 agitate in PBS for 1min.
- Remove excess buffer from sections.
- Drop commercial polymer anti-mouse HRP over the sections (enough to cover the tissue 50-100µl) incubate in a humid chamber following the manufacturer's instructions.
- Wash vigorously with PBS, agitate in PBS 1min.
- Remove excess buffer from sections and apply commercial Diaminobenzidine (DAB) following the manufacturer's instructions.
- Wash with water.
- Counter stain the sections briefly (2-3secs) with Harris hematoxylin diluted 50% in water.
- Wash with water, dehydrate and cover slip following normal laboratory protocols.

www.cenbimo.com
www.histosonda.com

Fabricante: CENBIMO S.L.
C/ Doctor Iglesias Otero, s/n
27004 Lugo, Spain